

► FiO₂ Sensors

12

MONTH
WARRANTY



MAX-1
R113P90

PRIMARY APPLICATIONS

- **Datascope** [Multinex Analyzer]
- **GE (Datex-Ohmeda)** [0237-2034-700 Excel Modulus, 5250 RGM, 5150 Handheld Monitor]
- **GE (Ohmeda)** [5120, 5100, 5250 RGM, 4700 Oxicap, 7800, 5125, 5150, Handheld Monitor, 0237-2034-700, Excel Modulus]
- **Hewlett Packard** [15201A]
- **Hudson** [5556]

SENSOR CROSS REFERENCE

- **Analytical Industries** [PSR-11-915]
- **Teledyne** [R15 #C-41131]

24

MONTH
WARRANTY



MAX-2
R107P20

PRIMARY APPLICATIONS

- **Drager** [6803290, 6850645, Narkomed, Fabius GS, Tiro]
- **Hewlett Packard** [15201A]
- **Invivo** [Magnitude]

SENSOR CROSS REFERENCE

- **Analytical Industries** [PSR-11-915-1]
- **Teledyne** [R31]

12

MONTH
WARRANTY



MAX-3
R109P11

PRIMARY APPLICATIONS

- **Mindray** [VT50 Ventilator]
- **Paragon** [Platinum SC430]
- **Penlon** [Prima SP2]

SENSOR CROSS REFERENCE

- **Analytical Industries** [PSR-11-60-03]
- **City Technology** [MOX-3]

12

MONTH
WARRANTY



MAX-7
R109P09

PRIMARY APPLICATIONS

- **Teledyne** [AX/MX300 Series,]

SENSOR CROSS REFERENCE

- **Analytical Industries** [PSR-11-77]
- **City Technology** [MOX-3]

► FiO₂ Sensors

14

MONTH
WARRANTY

MAX-9
R114P80

PRIMARY APPLICATIONS

- **HillRom** [6736140]
- **Hudson** [5500, 5590, 5577, 5584, 5595, 6477]

SENSOR CROSS REFERENCE

- **Analytical Industries** [PSR-11-55]
- **International Technologies** [M-13]
- **Ohio Medical** [655265]
- **Teledyne** [R13 #B-39636]

14

MONTH
WARRANTY

MAX-10
R112P10

PRIMARY APPLICATIONS

- **GE (Ohmeda)** [6050-0004-110, Modulus SE, 7900 Series, Aisys, Aestiva Series, SmartVent S-5, Aespire, Avance]

SENSOR CROSS REFERENCE

- **Analytical Industries** [PSR-11-915-4]
- **City Technology** [MOX-7]
- **EnviteC** [OOM-110]
- **International Technologies** [M-10]
- **Teledyne** [R30MED #B-729391]

12

MONTH
WARRANTY

MAX-11
R113P10

PRIMARY APPLICATIONS

- **Drager** [6803290, 6850645, Narkomed, Fabius GS, Evita Family, Babylog, Fabius Series]
- **Invivo** [Magnitude]
- **HillRom** [PSR-11-915-2]
- **Hudson** [5557]

SENSOR CROSS REFERENCE

- **Analytical Industries** [PSR-11-915-2]
- **City Technology** [MOX-16]
- **EnviteC** [OOM-201]
- **International Technologies** [M-02]
- **Teledyne** [R23 #C-611391]

12

MONTH
WARRANTY

MAX-11i
R113P06

PRIMARY APPLICATIONS

- **Drager** [Incubators, 18000 High Temp Application]

SENSOR CROSS REFERENCE

- **Philips (Respironics)** [NeoPAP]

► FiO₂ Sensors

24

MONTH
WARRANTY



MAX-12
R109P45

PRIMARY APPLICATIONS

- **Datascope** [0600-00-0002 (MAX-12F)]
- **Hamilton** [396008, 396009, Galileo, Raphael, Arabella, Aladdin, 51254, G5, Arabella]
- **Hudson** [5556]
- **IMT** [PF300]
- **Sensidyne** [Acutonic Hi Frequency Jet Ventilator, OOM102]
- **Spacelabs** [90513, Capnograph]
- **Taema** [YR049700, Cesar-2, Horus-4]

SENSOR CROSS REFERENCE

- **Analytical Industries** [PSR-11-917-M]
- **City Technology** [MOX-9]
- **EnviteC** [OOM-202]
- **International Technologies** [M-04]
- **Ohio Medical** [472062, 470644]
- **Siemens** [9004979-E347E]
- **Teledyne** [R24MED #B-62321]

14

MONTH
WARRANTY



MAX-12A
R109P48

PRIMARY APPLICATIONS

- **HillRom** [Isolette Incubator]
- **Philips (Respironics)** [NeoPAP]

SENSOR CROSS REFERENCE

- **Analytical Industries** [PSR-11-917-M1]
- **EnviteC** [OOM-106]

14

MONTH
WARRANTY



MAX-12C
R109P53

PRIMARY APPLICATIONS

- **Viasys** [Avea Ventilator, Infant SiPaP]

SENSOR CROSS REFERENCE

► FiO₂ Sensors

14

MONTH
WARRANTY

MAX-13
R115P10

PRIMARY APPLICATIONS

- **Covidien** [4-020933-00, 7820 Monitor, 7200 Vent]
- **Datascope** [Anestar S]
- **HillRom** [6735142]
- **Hudson** [5567]
- **Invivo** [9445 Sensor]
- **Philips (Respironics)** [Criterion Vent, Espirit Vent Lifecare (PLV-102)]
- **Philips-Respironics** [Criterion Vent, Espirit Vent Lifecare (#PLV-102), 37-079 Sensor]
- **Viasys** [66030 Sensor, 6400 Vent, 8400 Vent, VIP Vent]

SENSOR CROSS REFERENCE

- **Analytical Industries** [PSR-11-917-J]
- **City Technology** [MOX-2]
- **EnviteC** [OOM-102-1]
- **Ohio Medical** [406931, MiniOX 3000, MiniOX I, MiniOX III]
- **Teledyne** [R17MED #C-43690, R17MED #C-43690]

14

MONTH
WARRANTY

MAX-13-250
R125P07

PRIMARY APPLICATIONS

- **Covidien** [4-020933-00, 7820 Monitor, 7200 Vent]
- **Datascope** [Anestar S]
- **HillRom** [6735142]
- **Invivo** [9445 Sensor]
- **MiniOX I, MiniOX III** [406931, MiniOX 3000]
- **Philips (Respironics)** [Criterion Vent, Espirit Vent Lifecare (PLV-102)]
- **Philips-Respironics** [Criterion Vent, Espirit Vent Lifecare (#PLV-102), 37-079 Sensor]
- **Viasys** [66030 Sensor, 6400 Vent, 8400 Vent, VIP Vent]

SENSOR CROSS REFERENCE

- **Analytical Industries** [PSR-11-917-J]
- **Teledyne** [R17 #C-43690]

14

MONTH
WARRANTY

MAX-13C
R115P01

PRIMARY APPLICATIONS

- **BioMed Devices** [4433]

SENSOR CROSS REFERENCE

- **N/A**

► FiO₂ Sensors

14

MONTH
WARRANTY



MAX-14
R116P80

PRIMARY APPLICATIONS

- **Devilbiss** [O2ANA-101 Sensor]

SENSOR CROSS REFERENCE

- **Ohio Medical** [806572, MiniOX 1A]

16

MONTH
WARRANTY



MAX-16
R114P70

PRIMARY APPLICATIONS

- **Covidien** [4-072214-00, 840, 740, 760, PB2000 Tester]

SENSOR CROSS REFERENCE

- **Analytical Industries** [PSR-11-917M]
- **City Technology** [MOX-1]
- **EnviteC** [OOM-101-Air]
- **International Technologies** [M-15]
- **Ohio Medical** [716044, PB840]
- **Teledyne** [R-24 C-61340]

18

MONTH
WARRANTY



MAX-16HT
R114P73

PRIMARY APPLICATIONS

- **Versamed** [Vent]

SENSOR CROSS REFERENCE

- **N/A**

12

MONTH
WARRANTY



MAX-17
R116P10

PRIMARY APPLICATIONS

- **Hudson** [5558]
- **Penlon** [MTS]

SENSOR CROSS REFERENCE

- **Analytical Industries** [PSR-11-77]
- **City Technology** [MOX-4]
- **EnviteC** [OOM-105]
- **International Technologies** [M-07]
- **Ohio Medical** [804674]
- **Teledyne** [T-7]

► FiO₂ Sensors

12 MONTH WARRANTY



MAX-18
R116P40

PRIMARY APPLICATIONS

- **Hudson** [5569, 5525, 5528, 5529, 5550, 5552]

SENSOR CROSS REFERENCE

- **Analytical Industries** [PSR-11-33-1]
- **Ohio Medical** [485905]
- **Teledyne** [C1R #B-39022]

12 MONTH WARRANTY



MAX-19
R116P60

PRIMARY APPLICATIONS

- **Hamilton** [HM-01, Veolar, Amadeus, Arabella]
- **Hudson** [5565, 5568, C1, C-01]

SENSOR CROSS REFERENCE

- **Analytical Industries** [PSR-11-33-1]
- **EnviteC** [OOM-107]
- **International Technologies** [M-05]
- **Ohio Medical** [478841]
- **Teledyne** [T1 #A-37016]

12 MONTH WARRANTY



MAX-23
R116P06

PRIMARY APPLICATIONS

- **Teleflex** (Hudson) [5566]

SENSOR CROSS REFERENCE

- **Analytical Industries** [PSR-11-33-1]
- **City Technology** [MOX-1]
- **EnviteC** [OOM-102]
- **Teledyne** [R22 #C-44611]

16 MONTH WARRANTY



MAX-24
R109P15

PRIMARY APPLICATIONS

- **Stephan** [Stephanie Infant Ventilator]

SENSOR CROSS REFERENCE

- **N/A**

► FiO₂ Sensors

24

MONTH
LIFESPAN



PRIMARY APPLICATIONS

► N/A

SENSOR CROSS REFERENCE

► **Figaro** [KE-25]

MAX-25
R100P51-002

12

MONTH
LIFESPAN



PRIMARY APPLICATIONS

► **Hudson** [5811]

SENSOR CROSS REFERENCE

► N/A

MAX-27
R116P11

12

MONTH
LIFESPAN



PRIMARY APPLICATIONS

► **GE (Ohmeda)** [Giraffe]

SENSOR CROSS REFERENCE

► N/A

MAX-43
R109P07

12

MONTH
LIFESPAN



PRIMARY APPLICATIONS

- **Hudson** [5800 Monitor, 5801 Analyzer, 5803 Sensor]
- **Maxtec** [TecO₂ ME]

SENSOR CROSS REFERENCE

► N/A

MAX-47
R112P16

► FiO₂ Sensors

12

MONTH
WARRANTY

MAX-48
R112P18

PRIMARY APPLICATIONS

- **Analytical Industries** [All2001 Analyzer, All2002 Analyzer]

SENSOR CROSS REFERENCE

- **Analytical Industries** [All-11-60]

12

MONTH
WARRANTY

MAX-49
R116P13

PRIMARY APPLICATIONS

- **Hudson** [5804]

SENSOR CROSS REFERENCE

- N/A

12

MONTH
WARRANTY

MAX-50
R100P95-003

PRIMARY APPLICATIONS

- N/A

SENSOR CROSS REFERENCE

- **Figaro** [KE-50]

24

MONTH
WARRANTY

MAX-550E
R140P02

PRIMARY APPLICATIONS

- **Maxtec** [Maxblend 2, MaxO₂ ME]

SENSOR CROSS REFERENCE

- N/A

► FiO₂ Sensors

24
MONTH
WARRANTY



MAX-250
R125P01-002

PRIMARY APPLICATIONS

- **Fluke** [VT Mobile]
- **Maxtec** [OM-25A]

SENSOR CROSS REFERENCE

- **Analytical Industries** [PSR-11-75-KE1]
- **Teledyne** [R29]

24
MONTH
WARRANTY



MAX-250+
R125P02-011

PRIMARY APPLICATIONS

- **Maxtec** [MaxO₂+A]

SENSOR CROSS REFERENCE

- **Analytical Industries** [PSR-11-75-KE10+]

24
MONTH
WARRANTY



MAX-250A
R125P04-001

PRIMARY APPLICATIONS

- **Spacelabs**

SENSOR CROSS REFERENCE

- N/A

24
MONTH
WARRANTY



MAX-250E
R125P03-002

PRIMARY APPLICATIONS

- **CareFusion** [Sentry Blender]
- **Fluke** [VT Mobile]
- **Maxtec** [MaxO₂+AE, OM-25AE, OM-25ME, MaxBlend, MaxVenturi]
- **NewPort Medical** [e360 Vent, e500 Vent]
- **Viasys** [66045]

SENSOR CROSS REFERENCE

- **Analytical Industries** [PSR-11-917-J1]
- **Teledyne** [R29MED #C-70592]

► FiO₂ Sensors

12

MONTH
WARRANTY

MAX-250ESF
R115P85

PRIMARY APPLICATIONS

- **Maxtec** [MaxO₂+AE, OM-25AE, OM-25ME, MaxBlend, MaxVenturi]

SENSOR CROSS REFERENCE

- **Analytical Industries** [PSR-11-917-J1]
- **International Technologies** [M-25]

24

MONTH
WARRANTY

MAX-250TK
DC CONNECTOR
R125P11-001

PRIMARY APPLICATIONS

- **Viasys** [Vela Diamond Series, 16101]

SENSOR CROSS REFERENCE

- N/A

24

MONTH
WARRANTY

MAX-250K
PANDUIT CONNECTOR
R125P02-009

PRIMARY APPLICATIONS

- **Viasys** [Vela T2 (Older Models)]

SENSOR CROSS REFERENCE

- **Analytical Industries** [PSR-11-75-KE4]

24

MONTH
WARRANTY

MAX-250S
R125P16-001

PRIMARY APPLICATIONS

- Extra Stable Sensor for **MAX-250** Applications