



CAPSA HEALTHCARE MEDICAL CARTS



CAPSAHEALTHCARE

Avalo® Series Medical Carts



Avalo® Medical Carts

have long been a standard for promoting efficient, organized workflow, and providing reliable performance in facilities throughout the world.

Table of Contents

The Avalo Difference	3
Avalo Medical Cart Convenient Packages	4
Avalo Medical Specialty Carts	14
Avalo AC Anesthesia Cart	15
Product Specifications	16
Lock System Options	17
Configure Your Cart	18
M-Series Treatment Cart	20
Vintage Encore Treatment Cart	21



The Avalo® Difference

Avalo Medical Carts blend careful design, functionality, and logical organization to enhance user efficiency and support the delivery of superior patient care.

1. Durable, high-impact panels, built for long-term reliable performance. Stain, scratch and dent resistant.
2. Medications and supplies are secured with a selection of innovative locking systems including key, keyless or quick-access breakaway locking.
3. Flexible divider systems and removable drawer trays provide optimal visibility, grouping and organization of medications and supplies.
4. Seamless drawers, smooth panels and work surface promote simple cleaning.
5. Accessory options position equipment and supplies close at hand for easy access and improved efficiency, keeping the workspace clear.
6. Three drawer depths provide ample capacity to store and organize medications and supplies.
7. A lightweight ergonomic design and spread caster position provides stability control and ease of maneuverability.



Meeting the Needs for a Variety of Care Applications

Avalo Medical Carts can be deployed in a variety of care settings.

- Hospitals
- Assisted Living Facilities
- Surgery Centers
- Veterinary Clinics
- Physicians' Offices
- Dental Offices

Find the ideal Avalo Medical cart to fit your care application with new **easy-to-order convenient packages**, or **custom configure** a solution to meet specific workflow requirements.



Avalo® Series Medical Carts

Emergency Cart

The Avalo Series Emergency Cart blends careful design with functionality and logical organization to provide reliable performance and improved workflow in emergency departments. A unique breakaway locking handle offers optimal security and the swift, simple access essential to every emergency response. Stable maneuverability is built in and a full line of accessories enhance user efficiency.

Avalo Emergency Carts are available in convenient packages (shown below) or custom configurations.



Blue Emergency Cart
Standard Height
Shown with Accessory Package 2

Red Emergency Cart
Intermediate Height
Shown with Accessory Package 3




Quick-Access Lock unlocks with a simple and swift motion to quickly break the tamper-evident seal and gain access to all drawers.

Convenient Packages

Avalo Emergency Carts are offered in packages for your ordering convenience. In **two easy steps** choose a cart and an accessory package that best supports your workflow needs.

Step 1 Choose a Base Cart Option:

Choose a cart, then add an accessory package to configure the cart for your specific requirements.



Avalo® Series Medical Carts

Emergency Cart

Step 2 Choose a convenient accessory package:

Choose one of the following convenient accessory packages to support your facility's workflow, storage and organizational requirements. (Purchase in addition to cart)

Package 1

AM-EM-ACCPK1



Tamper Evident Seals (100)



Oxygen Tank Holder



Cardiac Board



I.V. Pole

Package 2

AM-EM-ACCPK2



Tamper Evident Seals (100)



Oxygen Tank Holder



Cardiac Board



I.V. Pole



Defibrillator Shelf



3" Drawer Dividers (2 Sets)

Package 3

AM-EM-ACCPK3



Tamper Evident Seals (100)




Oxygen Tank Holder



Cardiac Board



I.V. Pole



Defibrillator Shelf



3" Drawer Dividers (2 Sets)



Writing Surface



Sharps Container



Single Glove Box

Can't find what you're looking for?

If you require different drawer configurations or accessories not shown, please contact your Capsa Representative for help configuring a cart to meet your needs.



Avalo® Series Medical Carts

Pediatric Crash Cart

The **Avalo Series Pediatric Crash Cart** is designed for use with the proven clinical methodology for treating life-threatening emergencies in pediatric patients. When seconds matter, the color-coded drawers enable clinicians to quickly identify critical supplies when responding to emergencies.

The Avalo Pediatric Cart is available in convenient packages (shown below) or custom configurations.



Quick-Access Lock unlocks with a simple and swift motion to quickly break the tamper-evident seal and gain access to all drawers.

Convenient Packages

The Avalo Pediatric Crash Cart is offered in the standard height to support the proven clinical methodology for treating pediatric patients. In **two easy steps** choose the cart below and an accessory package to support your workflow needs.



Standard
AM-EM-STD-PED

Step 1 Choose the Standard Height Cart

Choose the standard cart, then add an accessory package to configure for your specific requirements.



Pediatric Crash Cart
Standard Height
Shown with Accessory Package 1



Avalo® Series Medical Carts

Pediatric Crash Cart

Step 2 Choose a convenient accessory package:

Choose one of the following convenient accessory packages to support your facility's workflow, storage and organizational requirements. (Purchase in addition to cart)

Package 1

AM-EM-ACCPK1



Tamper Evident
Seals (100)



Oxygen Tank Holder



Cardiac Board



I.V. Pole

Package 2

AM-EM-ACCPK2



Tamper Evident
Seals (100)



Oxygen Tank Holder



Cardiac Board



I.V. Pole



Defibrillator Shelf



3" Drawer Dividers
(2 Sets)

Package 3

AM-EM-PEDACCPK3



Tamper Evident
Seals (100)



Oxygen Tank Holder



Cardiac Board



I.V. Pole



Writing Surface



Defibrillator Shelf



3" Drawer
Dividers (8 Sets)

Can't find what you're looking for?

If you require different drawer configurations or accessories not shown, please contact your Capsa Representative for help configuring a cart to meet your needs.



Avalo® Series Medical Carts

Procedure / Treatment Cart

The Avalo Series Procedure/Treatment Cart offers ample storage capacity to accommodate the largest bulk items to the smallest medical supplies. The cart is adaptable to many different healthcare applications, including trauma, airway, dressing, cast, or general supply cart. A selection of proven lock systems provides secure storage.

Avalo Procedure/Treatment Carts are available in convenient packages (shown below) or custom configurations.



Seamless Drawer Trays are durable and prevent issues associated with pop rivets or sharp edges found in metal drawers.



*Procedure/Treatment Cart
Standard Height Shown with
Accessory Package 2*

Convenient Packages

The Avalo Procedure/Treatment Cart is offered in packages for your ordering convenience. In **two easy steps** choose the cart and accessory package to support your workflow needs.

Step 1 Choose a Base Cart Option:

Choose a cart, then add an accessory package to configure the cart for your specific requirements.



With Key Lock



Standard
AM-PT-STD-KEY



Intermediate
AM-PT-INT-KEY



Compact
AM-PT-CMP-KEY

Or



With Electronic Lock



Standard
AM-PT-STD-ELOCK



Intermediate
AM-PT-INT-ELOCK



Compact
AM-PT-CMP-ELOCK



Avalo® Series Medical Carts

Procedure / Treatment Cart

Step 2 Choose a convenient accessory package:

Choose one of the following convenient accessory packages to support your facility's workflow, storage and organizational requirements. (Purchase in addition to cart)

Package 1

AM-PT-ACCPK1



3" Drawer Divider
(1 Set)



Writing Surface



Waste with Lid

Package 2

AM-PT-ACCPK2



3" Drawer Divider
(1 Set)



Writing Surface



Waste with Lid



I.V. Pole



Sharps Container



Single Glove Box

Package 3

AM-PT-ACCPK3



3" Drawer Divider
(1 Set)



Writing Surface



Articulating Arm



Laptop Security Tray

Can't find what you're looking for?


If you require different drawer configurations or accessories not shown, please contact your Capsa Representative for help configuring a cart to meet your needs.



Avalo® Series Medical Carts

Anesthesia Cart

The Avalo Series Anesthesia Cart is designed to meet the unique needs of your anesthesia department. Flexible divider systems, removable drawer trays, and an array of accessories options promote orderly storage and organization. A choice of three innovative locking systems ensure that anesthesia medications and supplies are securely stored. Avalo Anesthesia Carts are available in convenient packages (shown below) or custom configurations.



Convenient Packages

The Avalo Anesthesia Cart is offered in packages for your ordering convenience. In **two easy steps** choose the cart and accessory package to support your workflow needs.

Step 1 Choose a Base Cart Option:

Choose a cart/lock option, then add an accessory package to configure the cart for your specific requirements.

Green Anesthesia Cart



Lock	Standard Height	Intermediate Height	Compact Height
Key Lock	AM-AN-STD-KEY-GR	AM-AN-INT-KEY-GR	AM-AN-CMP-KEY-GR
Elec Lock	AM-AN-STD-ELOK-G	AM-AN-INT-ELOK-G	AM-AN-CMP-ELOK-G
Auto Lock	AM-AN-STD-ALOK-G	AM-AN-INT-ALOK-G	AM-AN-CMP-ALOK-G

Blue Anesthesia Cart



Or

Lock	Standard Height	Intermediate Height
Elec Lock	AM-AN-STD-ELOK-B	AM-AN-INT-ELOK-B
Auto Lock	AM-AN-STD-ALOK-B	AM-AN-INT-ALOK-B




Avalo® Series Medical Carts

Anesthesia Cart


Step 2 Choose a convenient accessory package:

Choose one of the following convenient accessory packages to support your facility's workflow, storage and organizational requirements. (Purchase in addition to cart)

Package 1 AM-AN-ACCPK1



Package 2 AM-AN-ACCPK2



Package 3 AM-AN-ACCPK3



Can't find what you're looking for?
If you require different drawer configurations or accessories not shown, please contact your Capsa Representative for help configuring a cart to meet your needs.

Avalo® Series Medical Carts

Isolation Cart

The **Avalo Series Isolation Cart** provides ample drawer space to store and organize all of your isolation supplies, from gloves to gowns. Smooth panels, work surface and seamless drawers promote simple cleaning for infection prevention.

Drawers open fully for easy access to contents in the rear of the drawer and stay closed when the cart is in motion. Avalo's superior design offers durable performance for long-term value.

The Avalo Isolation Cart is available in convenient packages (shown below) or custom configurations.



Ample drawer space to store and organize all of your isolation supplies, from gloves to gowns.



*Isolation Cart
Standard Height*

Convenient Packages

The Avalo Isolation Cart is offered in packages for your ordering convenience.

In **two easy steps** choose the cart and accessory package to support your workflow needs.

Step 1 Choose a Base Cart Option:

Choose a cart, then add an accessory package to configure the cart for your specific requirements.



Standard
AM-IS-STD-NOLOK

Intermediate
AM-IS-INT-NOLOK

Compact
AM-IS-CMP-NOLOK

Avalo® Series Medical Carts

Isolation Cart

Step 2 Choose a convenient accessory package:

Choose one of the following convenient accessory packages to support your facility's workflow, storage and organizational requirements. (Purchase in addition to cart)

Package 1

AM-PT-ACCPK1



3" Drawer Divider
(1 Set)



Writing Surface



Waste with Lid

Package 2

AM-PT-ACCPK2



3" Drawer Divider
(1 Set)



Writing Surface



Waste with Lid



I.V. Pole



Sharps Container



Single Glove Box

Package 3

AM-PT-ACCPK3



3" Drawer Divider
(1 Set)



Writing Surface



Articulating Arm



Laptop Security Tray

Can't find what you're looking for?

If you require different drawer configurations or accessories not shown, please contact your Capsa Representative for help configuring a cart to meet your needs.



Avalo® Series Medical Carts

Specialty Carts

I.V. Therapy Cart

The Avalo Series I.V. Therapy Cart facilitates the organization of supplies for use in I.V. preparation. Keep everything in its place with optional drawer divider systems and storage bins. A large work surface offers ample room to prepare I.V. solutions. Expand your work area with a slide-out work surface, positionable on the right or left side of the cart.



Critical Care Cart

The Avalo Series Critical Care Cart is built to promote efficient workflow and task management in critical care environments. The durable frame and spread caster position ensure stability, while making the cart easy to move and control. A full complement of accessory options facilitates efficient workflow, and a choice of lock systems ensures the security of critical care supplies.



Medical/Surgical Cart

The Avalo Series Med/Surg Cart meets the demanding requirements of an ever-changing hospital environment. Every Med/Surg Cart accommodates the latest technology solutions and monitoring equipment to ensure efficient and streamlined workflow. Select from an array of drawer divider and tray solutions that offer flexibility for organizing supplies.



Avalo® Series Medical Carts

AC Anesthesia Self-Locking

The Avalo Series AC Anesthesia Cart offers the latest design and technology to meet the unique requirements of anesthesia storage, organization and security.

Flexible drawer configurations combined with optional drawer dividers and trays promote orderly storage.

An innovative keyless access and automatic relocking system provides secure and safe storage of anesthesia medications. Add another layer of security with an optional ID badge card reader, or optical bar code scanner and cart user audit.



Model ACM
Standard Height



Model ACS
Compact Height Shown
with optional accessories



Keyless Access with Automatic Relocking and optional security upgrades.



Three drawer depths provide ample capacity to store and organize medications and supplies.



Avalo® Series Medical Carts

Specifications

Standard Cart Models



Cart Specifications

	Surface Height	Depth	Width
Standard Height:	43.0"	24"	31"
	110 cm	61 cm	79 cm
Intermediate Height:	39.5"	24"	31"
	100 cm	61 cm	79 cm
Compact Height:	36.0"	24"	31"
	91 cm	61 cm	79 cm

Drawer Specifications

Height	Depth	Width
3"	17.5"	21.5"
7 cm	45 cm	55 cm
6"	17.5"	21.5"
16 cm	45 cm	55 cm
10"	17.5"	21.5"
25 cm	45 cm	55 cm

Specifications, availability & components are subject to change without notice.

AC Anesthesia Cart Models



Drawers can be configured to best meet your unique requirements.



AC Anesthesia Cart Dimensions

Base Width:	ACS:	25.5"	65cm
	ACM:	37.75"	96cm
Base Depth:	ACS:	24"	61cm
	ACM:	24"	61cm
Work Surface Height:	Compact:	39.5"	100cm
	Standard:	43"	110cm

Avalo® Medical Carts are Certified



Avalo® Series Security Features

Lock System Options

Lock System Choices

AutoLock™

Automatically secures a motionless unlocked medical cart with a user-defined relock time delay from 1-99 minutes. AutoLock permits up to 5,000 individual user codes to further enhance cart access tracking and security. SLA battery operated.

AutoLock security upgrades available.



Keyless Access with Manual Relocking

Simply enter a 4-digit access code, push down the lock lever and the cart is ready for access. To secure the cart, simply push in the lock lever. Clear visual indication of lock status is apparent from a distance.



Key Access with Manual Relocking

Featuring the Best® Core Removable Lock System, this basic yet proven lock system provides your medical cart advanced security and the flexibility to change keys easily.



Breakaway Lock

The Avalo Emergency Cart's unique breakaway locking handle is standard on every model, offers the security you require and ensures quick access to medications and supplies.



AutoLock™ Security Upgrades

Standard Host Software



Using Standard Host Software, authorized users can add, change or delete an access code for specific users.

Proximity Access



An advanced proximity card reader can be selected as an additional security access requirement.

Keyless Narcotics Storage



The keyless narcotics storage feature augments security to further optimize control and accountability.

Available for Avalo AC Cart only

Custom Configured Carts

Build Your Custom Cart in 4 Steps

If specific drawer configurations or accessories are required to meet unique care applications, Avalo carts can be fully configured in a series of simple steps. Please contact your Capsa Representative for help configuring a cart to meet your needs.

Step 1 Select your cart height & drawer configuration



Choose drawer configuration

Choose a drawer configuration combined with optional drawer dividers and trays to meet your storage needs.

Step 2 Select your lock system (see page 17 for descriptions)



Step 3 Select your color and drawer I.D. label

Choose from our palette of colors



Choose from a selection of drawer ID labels

- Emergency
- Anesthesia
- Med/Surg
- Treatment
- Isolation
- Procedure
- IV Therapy
- Emergency Preparedness

Step 4 Select your accessories (see page 19)

Customize your cart with optional accessories to support the unique storage, organization and efficient workflow required.

Avalo® Series Medical Carts

Options & Accessories



Contact Capsa Healthcare to inquire about additional Accessories & Options



M-Series Medical Carts

Treatment Cart

The M3-TC Treatment Cart offers a lightweight aluminum frame, with large storage capacity and easy maneuverability. Smart design provides plenty of room for supplies, and drawers are customizable for optimal organization. The M3-TC improves workflow and medical treatment preparation.



Treatment Cart Features

- Large capacity in a combination of shallow, medium, and deep drawers
- Integrated top handle and full-width drawer handles
- Reliable and easily changeable lock mechanism
- Durable design for long-term value
- 10-year limited warranty

Accessories & Options



Waste w/ Lid



Organizer Tray



Drawer Flex Dividers



Chair Rail Bumpers



Writing Surface



Sharps Container



Externals Drawer Dividers



Arm Mount And Laptop Tray

Contact Capsa to inquire about additional Accessories & Options

Cart Specifications

Cart Base 39.4"h x 26.75"d x 29.5"w

Drawer Specifications

	Height	Depth	Width
Four Inch	4.25" / 11 cm	18.5" / 47 cm	22" / 56 cm
Six Inch	6.18" / 16 cm	18.5" / 47 cm	22" / 56 cm
Nine Inch	9.18" / 23 cm	18.5" / 47 cm	22" / 56 cm

M-Series Accent Colors



Blue-Gray



Platinum



Teal



Cranberry



Reflex Blue



Encore Series Medical Carts

Treatment Cart

The Vintage Encore Treatment Cart combines furniture styling with lightweight aluminum construction and impressive storage capacity. Smart design provides plenty of room for supplies, and drawers are customizable for optimal organization. The Vintage Encore Treatment Cart is ideal for modern decor, available in four finishes.



Treatment Cart Features

- Elegant styling to suit any facility's décor, in four furniture-style finishes
- Large capacity in a combination of shallow, medium, and deep soft-close drawers
- Simple-to-push casters
- No-bounce extended work surface
- Durable design for long-term value
- 10 year limited warranty

Accessories & Options



Waste w/ Lid



Organizer Tray



Arm Mount And Laptop Tray



Glove Box



Writing Surface



Sharps Container



10" Narcotic Box



Supply Dividers



Externals Drawer Dividers

Contact Capsa to inquire about additional Accessories & Options

Cart Specifications

Cart Base 40.1"h x 25.2"d x 28"w

Drawer Specifications

Drawer	W x D x H (inch)	W x D x H (cm)
Four Inch	23 x 18.5 x 4.25	58 x 47 x 10.8
Six Inch	23 x 18.5 x 6	58 x 47 x 15.2
Nine Inch	23 x 18.5 x 9	58 x 47 x 22.9

Available Finishes



White Ash



Boston Cherry



Antique Walnut



Chatsworth Oak



Additional Information

Contact Capsa Healthcare for more information about our wide selection of Medical Cart solutions that promote orderly storage and organization where it matters most. Adaptable to fit any application, Capsa's procedural and supply solutions are built to meet the unique needs of facilities worldwide.



Website

Visit the Capsa website to learn more about Capsa's wide selection of procedural and supply solutions.

www.CapsaHealthcare.com



Email

Email us to request more information about Capsa products or to inquire about regional representation or distribution partners in your area.

Info@CapsaHealthcare.com



Phone

If you have specific questions and wish to speak with someone directly, contact Capsa Healthcare and ask for Customer Service.

800.437.6633 International: +1.614.864.9966

Look for us on:



Certified ISO
27001:2013



Certified ISO
9001:2015



CapsaHealthcare.com | 800.437.6633

8170 Dove Parkway
Canal Winchester, Ohio 43110

Carro de de dosis unitaria para 18 pacientes



Seguro Digital



Residuo Interno



Extra Swivel W/ Brake

Tracking Caster